

Counting Protoplast with the LUNA-FX7™ Automated Cell Counter: Optimal Fluorescent Dye Combinations

Introduction

Protoplasts are plant cells with removed cell walls through enzymatic or mechanical removal of the cell wall. These spherical, cell-wall-free plant cells exhibit unique characteristics, including totipotency, making them versatile tools in plant science. They serve as an experimental system that allows to explore the structure, chemistry, and function of plant cells.

Traditionally, viability assessment involved culturing protoplasts until they developed into complete plants, but this method could not provide immediate viability assessment. The LUNA-FX7™ Automated Cell Counter now enables rapid viability determination through double staining with distinct-colored fluorescent dyes. However, it is essential to understand how the LUNA-FX7™ Automated Cell Counter performs with different dyes and to adjust analysis parameters accordingly. This study aims to identify optimal dye combinations and parameters for assessing protoplast viability using the LUNA-FX7™ Automated Cell Counter.

Protoplast Isolation and Staining

Protoplast Isolation

1. Prepare 5 g of seedlings and wash them once with ethanol and twice with distilled water (DW).
2. Gently pat dry the seedlings with a paper towel and place them in a Petri dish. Use scissors to finely chop them.
3. Add 25 ml of digestion buffer to the finely chopped seedlings, mix well, and transfer to a 50 ml tube. Cover with foil to protect from light.
4. Incubate on a rotator at 20 rpm for 6 hours.
5. After incubation, add 20 ml of wash buffer to the tube, gently mix, and filter through a 100 µm strainer.
6. Centrifuge at 100 g for 3 minutes.
7. Remove the supernatant, add 20 ml of wash buffer to the tube, gently mix, and filter through a 40 µm strainer.
8. Centrifuge at 100 g for 3 minutes.
9. Remove the supernatant and resuspend the pellet in 3 ml of wash buffer.

Fluorescence Staining

1. Mix:
 - 18 μ L of protoplast cells
 - 2 μ L of dyes prepared by mixing green and red fluorescent dyes
2. Load 10 μ L of stained cells.
3. Perform analysis using the LUNA-FX7™.

*Note: Protoplasts have diverse sizes and may exhibit varying fluorescent intensities. Please make adjustments as needed.

Table 1. The suggested parameter settings for protoplast counting of the LUNA-FX7™ on Fluorescence Cell Counting mode

Basic		Advanced	
Counting Mode	Fluorescence Cell Counting	Counting Mode	Fluorescence Cell Counting
GF exposure level	5	GF exposure level	5
RF exposure level	9	RF exposure level	9
Cell size calculation	BF	Min. cell size	3 μ m
Min. cell size	3 μ m	Max. cell size	70 μ m
Max. cell size	70 μ m	Declumping sensitivity	5
GF threshold level	5	Min. FL intensity	0
RF threshold level	5	Min. roundness	3
Dilution factor	1.11	Dilution factor	1.11

Fluorescent Dyes for Viability Measurement

We have selected five commonly used fluorescent dyes for viability assessment:

Dye	Properties	Colors
Calcein AM	Membrane-permeable Esterase-activated dyes	Green
Fluorescein diacetate (FDA)		
Acridine Orange (AO)	Membrane-permeable nuclear dye	Red
Propidium Iodide (PI)	Membrane-impermeable nuclear dye	
Ethidium Homodimer-1 (EthD-1)		

In principle, live cells should emit a green fluorescence when stained with membrane-permeable dyes while dead cells are expected to produce a red signal when stained with membrane-impermeable dyes. These dyes function through distinct mechanisms. For instance, nuclear dyes like AO, PI, and EthD-1 become significantly brighter when bound to DNA. Additionally, both Calcein AM and FDA require esterase activity to generate a green fluorescent signal which serves as an indicator of cell vitality for assessing cell health. In our experiments, we paired these green and red dyes to determine the viability of protoplasts, thus providing a comprehensive assessment of cell viability.

Optimal Fluorescent Dye Combinations for Protoplast Viability Assessment

FDA/PI or FDA/EthD-1 is the most effective choice for protoplast viability assessment (Figure 1). Both PI and EthD-1 exhibited the ability to stain protoplast nuclei; however, one may need to adjust the RF exposure level to 9 to detect sufficient signals (Table 1).

Among the green dyes we tested, only FDA yielded a bright and reliable signal with no need of adjusting GF exposure level. Both FDA and Calcein AM rely on esterase activity to produce signals, but Calcein AM exhibited almost no signal. However, it's important to be aware that FDA can generate high background noise during prolonged incubation periods. We recommend performing cell counting immediately after staining to achieve optimal results in viability assessment. If the signal is excessively high, adjusting the GF exposure level can be applied. Moreover, we found a significant decrease in AO intensity although AO is generally expected to stain all cells.




Figure 1. The staining results of isolated protoplast. Protoplasts were effectively stained with FDA/PI or FDA/EthD-1, whereas AO/PI or EthD-1 and Calcein AM/PI or EthD-1 exhibited no green signals.

FL: Images merged from both green and red fluorescence channels.

Tag: Composite images of all channels, fluorescent and brightfield, with identified objects marked using red and green circles. Red circles indicate dead cells, while green circles represent live cells.

Buffer-Dependent Performance of Green Fluorescent Dyes

Despite of AO's known ability to stain all cells regardless of viability, the unexpected low AO signal when used with protoplasts led to further investigation. We conducted experiments using green dyes on mammalian cells, including U937 cells, which are known to be stained by AO. The buffer was changed to either PBS or a wash medium to assess the impact on dye signals, and the default protocol was used for analysis. Indeed, there was a significant reduction in AO signal intensity in the wash medium compared to PBS (Figure 2A). Moreover, Calcein AM exhibited more consistent signals when the washing buffer was used rather than PBS, and FDA showed a minor decrease in overall intensity when the washing buffer was used (Figure 2B and 2C).

Although the precise factors contributing to the differences in dye performance remain unclear, the differences in osmolarity and pH between the wash medium and PBS may influence cells and dye performance. For example, it's important to note that AO can display varying colors at different pH levels. The difference in osmolarity made a reduction in cell size approximately 13 μm to 1-2 μm smaller. Considering these factors, these dyes can be influenced by diverse mechanisms based on the chemical conditions surrounding the cells.




Figure 2. The impact of different buffers on the performance of AO (A), Calcein AM (B), and FDA (C) on U937 cells. The AO signal decreased in wash medium, while Calcein AM exhibited improved performance in wash medium. The performance of FDA was not influenced by the type of buffer.

BF: Images captured in brightfield.

FL-Tag: Composite fluorescent channel images with identified objects marked using red and green circles. Red circles indicate dead cells, while green circles represent live cells.

Conclusion

FDA/PI or FDA/EthD-1 dyes are the most effective dyes that effectively stain protoplasts. While FDA consistently produces a reliable signal, it is important to conduct cell counting immediately to minimize potential background noise during incubation. Calcein AM showed discrepancies in performance depending on the cell types and buffers. Interestingly, the AO signal exhibits a significant reduction when the wash medium is used. Additionally, both PI and EthD-1 are effective in staining protoplasts, but it requires an RF exposure level adjustment from 5 to 9. In summary, choosing the appropriate dye combinations and optimal analysis protocols of the LUNA-FX7™ Automated Cell Counter serve as a great method for protoplast viability assessment.